

O and P¹ • 1.800.408.3598

Friddle's • 1.800.528.9339

Spyder Strong!

SPYDER
TECHNOLOGIES

BRAND-NEW, EXCLUSIVE AND
Spyder Strong!

O and P¹, Waterloo, Iowa, and Friddle's, Honea Path, S.C., are the exclusive manufacturers of the new Spyder Technologies custom socket.

Up to 60 percent lighter, and stronger than any manufactured socket today, the braided socket is being used by some of the world's elite amputee athletes and soldiers.

Carbon fibers are braided individually to the patient's shape - not simply laminated with manufactured carbon stockinet - and are manufactured to meet stringent ISO standards. The technician has the ability to control the braid while it is being woven directly on the positive plaster mold or foam carving.

Spyder is

- Clinically engineered
- Engineer tested for strength
- Manufactured to exceed ISO standards
- Comparable in price to similar prostheses
- A superior product that has been three years in development
- Can make a 1+ pound difference in a prosthetic leg

